

## Composition

### P07-03100 ITS solution I (100x)

|                  | Components                                     | mg/L     |
|------------------|------------------------------------------------|----------|
| Inorganic Salts  | Calcium chloride anhydrous                     | 200.00   |
|                  | Magnesium sulfate x 7H <sub>2</sub> O          | 200.00   |
|                  | Potassium chloride                             | 400.00   |
|                  | Sodium chloride                                | 6,800.00 |
|                  | Sodium dihydrogen phosphate x H <sub>2</sub> O | 140.00   |
|                  | Sodium selenite                                | 0.500    |
| Other Components | D(+)-Glucose anhydrous                         | 1,000.00 |
|                  | Insulin human rec.                             | 500.00   |
|                  | apo-Transferrin human                          | 500.00   |