

Composition

P04-02504 **SILAC-RPMI 1640**
w/o: L-Glutamine
w/o: L-Arginine
w/o: L-Lysine
w: 2.0 g/L NaHCO₃

	Components	mg/L
Inorganic Salts	Calcium chloride x 2H ₂ O	62.30
	Potassium chloride	400.00
	Magnesium sulfate anhydrous	48.83
	Sodium chloride	6,000.00
	di-Sodium hydrogen phosphate	800.00
	Sodium nitrate	72.00
Other Components	D(+)-Glucose anhydrous	2,000.00
	Glutathion (red.)	1.00
	Phenol red	5.00
Amino Acids	L-Arginine x HCl	0.00
	L-Asparagine x H ₂ O	50.00
	L-Aspartic acid	20.00
	L-Cystine	50.00
	L-Glutamine	0.00
	L-Glutamic acid	20.00
	Glycine	10.00
	L-Histidine free base	15.00
	L-Hydroxyproline	20.00
	L-Isoleucine	50.00
	L-Leucine	50.00
	L-Lysine x HCl	0.00

	L-Methionine	15.00
	L-Phenylalanine	15.00
	L-Proline	20.00
	L-Serine	30.00
	L-Threonine	20.00
	L-Tryptophan	5.00
	L-Tyrosine	20.00
	L-Valine	20.00
Vitamins	p-Aminobenzoic acid	1.00
	D-(+)-Biotin	0.20
	D-Calcium pantothenate	0.25
	Choline chloride	3.00
	Folic acid	1.00
	myo-Inositol	35.00
	Nicotinamide	1.00
	Pyridoxine x HCl	1.00
	Riboflavin	0.20
	Thiamine x HCl	1.00
	Vitamine B12	0.005
	NaHCO ₃	2,000.00