

**Datasheet****Bovine Serum Albumin****30% Solution**

Product	Description	Catalogue-No.	Size
<b>BSA 30% Solution</b>	<b>Bovine serum albumin, Fraction V, 30% Solution</b>	<b>P06-1403020</b>	<b>20 ml</b>
		<b>P06-1403100</b>	<b>100 ml</b>
		<b>P06-1403500</b>	<b>500 ml</b>

**Product description**

BSA 30% solution is a solution prepared from BSA, Fraction V. The molecular weight of BSA is 66,430 Da. BSA is a single polypeptide chain consisting of about 583 amino acid residues and no carbohydrates. At pH 5-7 it contains 17 intrachain disulfide bridges and 1 sulfhydryl group.

Albumins serve as additive proteins for tissue cultures. They are the main protein component in serum and are added to cell culture media to increase the stability of cell membranes and to bind possibly toxic elements. It also binds fatty acids, bilirubin and hormones.

**Storage conditions**

Storage: -20 °C

Stability: 2 years from date of production

Filling: 20 ml, 100 ml, 500 ml, other sizes on request

**Composition**

BSA 30% solution is a non-preserved solution produced from Bovine Serum Albumin, Fraction V.

**Suitability**

BSA 30% solution high monomer is suitable for use as a protein stabilizer for antibody and enzyme storage.

**Special advantages**

- Highly purified solution absent of chemical and microbial contamination
- Bovine serum albumins increase the stability of cell membranes and bind possibly toxic elements
- Ready to use solution
- Manufactured under conditions certified to ISO 9001

**Technical support**

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email ([info@pan-biotech.com](mailto:info@pan-biotech.com)) or phone +49-8543-601630.

---

FOR RESEARCH USE ONLY! Not approved for human or animal diagnostic or therapeutic procedures.