

**Datasheet****Sterile Water, for Cell Culture****Media Supplement**

Product	Description	Catalogue-No.	Size
Sterile Water	Water for Cell Culture	P04-99025	2.5 L
		P04-991000	1000 ml
		P04-991500	500 ml
		P04-99200B	20 L
		P04-995000	5 L

**Product description**

PAN-Biotech's Sterile Water for Cell Culture is purified water which is prepared by ultrafiltration, reverse osmosis, electrodeionization and sterile filtration.

**Storage conditions and stability**

Storage: +2 °C to RT

Stability: 3 years from date of production

Size: 500 ml, 1000 ml, 2.5 L, 5 L, 20 L, other sizes on request

**Release Criteria**

Endotoxin < 1 EU/ml

Sterility pass

**Suitability**

PAN Biotech Sterile Water for Cell Culture is intended for research use, not for human or animal diagnostic or therapeutic procedures including e.g. parenteral administration.

**Special Advantages**

- Highest quality standards
- Offered in flexible packaging solutions
- Lot-specific Certificate of Analysis (CoA) for each lot of water produced

**Technical support**

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email ([info@pan-biotech.com](mailto:info@pan-biotech.com)) or phone +49-8543-601630.

---

FOR RESEARCH USE ONLY! Not approved for human or animal diagnostic or therapeutic procedures.