

## Composition

**P07-03400 ITS solution IV (100x)**  
**P07-03410**

|                  | Components                                     | mg/L        |
|------------------|------------------------------------------------|-------------|
| Inorganic Salts  | Potassium chloride                             | 400.00      |
|                  | Sodium chloride                                | 6,800.00    |
|                  | Sodium dihydrogen phosphate x H <sub>2</sub> O | 140.00      |
|                  | Sodium selenite                                | 0.500       |
| Other Components | BSA                                            | 50,000.00   |
|                  | D(+)-Glucose anhydrous                         | 1,000.00    |
|                  | Linoleic acid-Oleic Acid-Albumin (100x)        | 5.00 (ml/L) |
|                  | Insulin human rec.                             | 1,000.00    |
|                  | apo-Transferrin human                          | 550.00      |