

Datasheet

William's Medium E**Cell Culture Medium**

Product	Description	Catalogue-No.	Size
William's Medium E	Cell Culture Medium	P04-29510	500 ml

Product description

The William's Medium E is used for long-term cultivation of adult rat liver epithelial cells. It can also be used for cultivating other mammalian cells and is usually supplemented with 5-10 % FBS to achieve long term cell cultivation.

Storage conditions

Storage: 2-8°C

Stability: 2 years from date of production

Composition

William's Medium E

w/o: L-Glutamine

w/o: Phenol red

w: 2.24 g/L NaHCO₃

Detailed composition upon request.

Performance Testing

William's Medium E is tested for pH, osmolality and sterility. A Certificate of Analysis is available upon request.

Technical support

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email (info@pan-biotech.com) or phone +49-8543-601630.

FOR RESEARCH USE ONLY! Not approved for human or animal diagnostic or therapeutic procedures.